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**steriCLIN®**

## Hersteller

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**Product**

See-through packaging made of Tyvek® and film

**Article Number**

3FKFB250xxx, 3FKFS250xxx

**Intended Purpose**

Pre-manufactured sterile barrier system to be used in packaging systems for terminally sterilised medical devices.

Purpose of these see-through packaging is to enable sterilisation, physical protection, maintenance of sterility and aseptic presentation of packed medical devices.

**Medical Device**

**According to the intended purpose, the product is a medical device.**

Conformity for medical devices of class I is declared according to Medical Device Regulation (EU) 2017/745.

The see-through packaging is CE-marked.

**Intended Users**

Intended users is educated and trained staff working in the Central Sterile Supply Department of health care facilities.

**Clinical Benefit**

These sterile barrier systems are not directly applied to patients and do not have a medical indication. Prevention of infections is the clinical benefit.

**Performance Characteristics**

The performance characteristics do result from compliance with the product related standards for sterile barrier systems. The products comply with applicable standards ISO 11607 and EN 868-5. Imprinted chemical indicators comply with standard ISO 11140-1 (Type 1).

**Intended Method of Sterilisation**

The see-through packaging made of Tyvek® and film is supplied in non-sterile condition. In application, they are suitable for the methods of sterilisation as follows:

**VH2O2 | IRRAD | EO | FORM**

**Important Instruction**

Serious incidents with this see-through packaging shall immediately be reported to the manufacturer and to the competent national authority relevant for medical devices.

**Warnings**

- **Do not re-use!** This see-through packaging is intended for single use. Sterile barrier functionality may be compromised in case of multiple use.
- Sterile barrier functionality may be compromised by influence of damaging substances and chemicals.
- For sterilisation method other than intended, the proper functionality of sterile barrier cannot be guaranteed.
- If items with sharp edges and tips are packed, damages of the sterile barrier system may happen and result in a loss of sterile barrier.

**Precautions**

**Do not use when the sterile barrier system or its packaging is damaged!** In general, contact to chemicals, detergents or disinfectants shall be avoided. Before sealing, the sterile barrier systems shall be conditioned for 24 hours under prevailing climatic conditions.

**Safe Disposal**


Non-contaminated sterile barrier systems can be disposed of together with household waste.

# Instructions for Use


## Storage Conditions and Stability


Recommended storing condition is at temperatures between 15°C and 30°C and at relative humidity level between 35% and 60%. **Keep dry! Keep away from sunlight and heat!** The shelf life of see-through packaging until the use for packaging of a medical device, sealing and sterilisation does cover 5 years. Subsequently, the stability of the sterile barrier is guaranteed for further 5 years of storing. Other limitations of storage period (e.g. recommendations of DIN 58953-8) have to be noted with respect to the storage conditions.


## Use

Before packaging, tips and sharp edges of products should be covered in a suitable manner.

During packaging of products, the later opening direction (peel direction) shall be considered.

The use of the filling area shall not exceed 75% of available filling space.


A distance of > 3 cm shall be left between the cutting and the seal.

When sealing on common rotary sealers, a recommended sealing temperature from 115°C to 135°C shall be used.

Tyvek® shall always contact to Tyvek® when producing double packaging.


Markings shall be placed on the film side only, and always outside of the filling area.

To read the chemical indicators, the colour shall be compared with information as printed on the packaging.

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