

Bowie & Dick single use test pack - ISO 11140-4

The utmost challenge for your steriliser

Correct functionality of a steriliser is key when it comes to ensuring excellent working results in the central sterile supply department (CSSD). The stericlin® Bowie & Dick single use test pack enables daily function checks and quality control of sterilisers.

Complete ease of handling: Place the test pack in your steriliser, run the Bowie & Dick test programme, remove the test sheet, check the result, and that's it. The indicator on the included test sheet changes from blue to pink, and contains no heavy metals. After use, the colour of the indicator no longer changes, ensuring total reliability of documentation.

The stericlin® single use test pack for performing daily Bowie & Dick tests for steam sterilisers was awarded the prestigious Kitemark quality certification of the British Standards Institute (BSI). This confirms the product's compliance with the stringent DIN EN ISO 11140-4 standard requirements by an internationally accepted, independent testing institute.

Best quality in compliance with standards and fulfilment of ISO 11140-4 requirements can and must be expected from a test tool like the stericlin® Bowie & Dick single use test pack. This is only possible when using a specially equipped steriliser as implemented by the development laboratory at VP. BSI testers verified the fulfilment of prerequisites for this test as well as correct execution, on site and in person.

The Kitemark of the BSI is one of the oldest symbols of quality ever and is especially awarded to products and services where safety plays a crucial role, subject to the most stringent conditions.



Bowie & Dick single use test pack

Article-No.	Description	Pack (pcs.)
3FSKS611102	Bowie & Dick single use test pack	40